

Oxivir Spray E2n

Cleaner & disinfectant, oxygen-based

E2n

Description

Diversey Oxivir Spray is a ready-to-use oxygen based liquid cleaner disinfectant for the cleaning and disinfection of all water-resistant hard surfaces.

Key properties

Hydrogen peroxide based acidic formulation
Effective against a wide spectrum of micro-organisms
Effective cleaner
Fragrance-free
Based on patented technology; accelerated hydrogen peroxide

Benefits

Safe to use, efficiently hinders lime scale build-up
Passes new European norms: Virucidal, Fungicidal and Bactericidal
Cleans and disinfects in one step, reduction of working time
Streak-free cleaning
Suitable for use in all kind of application areas, eg in high risk areas in Healthcare
Environmentally more friendly than common disinfectants

Use instructions

Dosing:
Use undiluted

Application:

Spray cleaning and disinfection:

1. Remove loose dirt (for heavy soiled areas, always pre-clean with a JohnsonDiversey cleaner and rinse with water)
2. Spray onto surfaces
3. Clean with clean cloth, pre-wetted with product
4. Spray again and leave surface wet for at least 5 minutes.
5. Rinse surfaces thoroughly with clean water and allow to air dry

Important:

Do not mix with other products. Do not use on water- or acid sensitive surfaces.



Technical data

Appearance:	clear, colourless liquid
Relative density [20°C]:	1.0
pH-value in use:	1.5

The above data is typical of normal production and should not be taken as a specification

Safe handling and storage information

Full guidance on the handling and disposal of this product is provided in a separate Material Safety Data Sheet. Store in original container. Avoid extremes of temperature. Reacts with oxydisable substances.



Product compatibility

Under recommended conditions of use **Diversey Oxivir Spray** is suitable for use on most materials commonly encountered.

Microbiological data

Diversey Oxivir Spray cleaner disinfectant passed EN 1276 at a contact time of 5 minutes at 20°C under dirty conditions and in hard water.

Diversey Oxivir Spray cleaner disinfectant passed EN 13727 at 20°C at a contact time of 5 minutes

Diversey Oxivir Spray cleaner disinfectant tested according to EN 14476 showed to be effective against Poliovirus Type 1 and Adenovirus Type 5 at a contact time of 5 minutes under clean conditions, without interfering substances.

Available Pack size

6 x 750ML

Diversey Eastern and Central Africa Ltd

Kaptagat Road, Loresho

P.O Box 41939 - 00100 Nairobi

Tel. +254-20-4224000 3008000

Fax:+254(20)4224888 / 4180481

www.diversey.com